





A Family of DC Magnetic Sensors to Best Meet Your Requirements

medSAFE™

Sensor options:

<p>Model 90 6DOF Sensor</p>		<p>Model 90 6DOF Sensor</p> <ul style="list-style-type: none"> • Sensor OD = 0.9 mm • Sensor Length = 7.25 mm • Cable OD = 0.6 mm • Cable length = 6.6 ft. (2.0 m)
<p>Model 130 6DOF Sensor</p>		<p>Model 130 6DOF Sensor</p> <ul style="list-style-type: none"> • Sensor OD = 1.5 mm • Sensor Length = 7.6 mm • Cable OD = 1.2 mm • Cable length = 6.6 ft. (2.0 m)
<p>Model 180 6DOF Sensor</p>		<p>Model 180 6DOF Sensor</p> <ul style="list-style-type: none"> • Sensor OD = 2 mm • Sensor Length = 9.9 mm • Cable OD = 1.2 mm • Cable length = 6.6 ft. (2.0 m)
<p>Model 800 6DOF Sensor</p>		<p>Model 800 6DOF Sensor</p> <ul style="list-style-type: none"> • Sensor OD = 7.9 mm • Sensor Length = 19.8 mm • Cable OD = 3.8 mm • Cable length = 6.6 ft. (2.0 m)

