





**A Family of DC Magnetic Sensors to Best Meet Your Requirements**

# medSAFE™

## Sensor options:

<p><b>Model 90 6DOF Sensor</b></p>		<p><b>Model 90 6DOF Sensor</b></p> <ul style="list-style-type: none"> <li>• Sensor OD = 0.9 mm</li> <li>• Sensor Length = 7.25 mm</li> <li>• Cable OD = 0.6 mm</li> <li>• Cable length = 6.6 ft. (2.0 m)</li> </ul>
<p><b>Model 130 6DOF Sensor</b></p>		<p><b>Model 130 6DOF Sensor</b></p> <ul style="list-style-type: none"> <li>• Sensor OD = 1.5 mm</li> <li>• Sensor Length = 7.6 mm</li> <li>• Cable OD = 1.2 mm</li> <li>• Cable length = 6.6 ft. (2.0 m)</li> </ul>
<p><b>Model 180 6DOF Sensor</b></p>		<p><b>Model 180 6DOF Sensor</b></p> <ul style="list-style-type: none"> <li>• Sensor OD = 2 mm</li> <li>• Sensor Length = 9.9 mm</li> <li>• Cable OD = 1.2 mm</li> <li>• Cable length = 6.6 ft. (2.0 m)</li> </ul>
<p><b>Model 800 6DOF Sensor</b></p>		<p><b>Model 800 6DOF Sensor</b></p> <ul style="list-style-type: none"> <li>• Sensor OD = 7.9 mm</li> <li>• Sensor Length = 19.8 mm</li> <li>• Cable OD = 3.8 mm</li> <li>• Cable length = 6.6 ft. (2.0 m)</li> </ul>

